

USE OF THE

TK-4

COMPRESSION STRAP TOURNIQUET



Tear open the package by using one of the indicated tear notches and remove the product. Grab the product by one of the S-hooks. Then follow the instructions below.

1

Place the TK-4 by the hook end over the injured limb above the wound area. Lay the strap over the limb with the hook to the outside of the limb.

2

Wrap the elastic strap under the limb and hook it through the S-hook on the limb. Make sure the elastic is secure in the hook. Begin pulling on the strap and keep a constant tension on the elastic.

3

Holding the TK-4 against the limb, pull and tighten the strap down to secure the TK-4 and begin compression against the limb. Hold tight to the strap and begin winding it around the limb, continuing to pull the strap tight as you wrap. **DO NOT** let go of the strap while winding; this will result in the hook and strap unraveling and may cause injury.

4

Once the tourniquet has been wrapped around the limb, tuck the remaining S-hook under the elastic wrap to secure the tourniquet. If you need to adjust the tourniquet, pull out the S-hook and unwrap the strap.



Frequently check any style tourniquet for continued occlusion as loosening can occur after application.



THIS PRODUCT IS SINGLE USE ONLY. DO NOT REUSE.

This product is a medical device that may require specific training regarding proper use. Please follow any local guidance or governing authority regarding the use of this device.

