



Wound Seals

H&H Med Corp Continues to Lead in the Innovation of Pneumothorax and Hemothorax Care.



H*VENT

Fail-Safe Design

6 ports allow for redundancy

Multi-directional Design

Allow the fluids to “drain with gravity”

Large Pull Tab

Help with application and repositioning

SINGLE PACK (HVS01)

NSN: 6510-01-671-7225

Dimensions: 7.75" x 4.5" x 0.25"

Weight: 1.4 oz.

TWIN PACK (HVS02)

Dimensions: 9.25" x 8.5" x 0.25"

Weight: 2.9 oz.

WOUND SEAL

Rounded Edges

Reduce friction with clothing

Long Pull Tab

Allow venting, help with repositioning and removal

Large Coverage Area

Allow hydro-gel adhesive to stick through blood/debris

HHWSK02

NSN: 6510-01-573-0300

Size: 8.25" x 6.25" x 0.25"

Weight: 2.5 oz.

Sold in the H&H Wound Seal kit:

Individual kits or in cases of 25

Tested to MIL-STD-810G

Shelf Life: 5 years

BOLIN CHEST SEAL

Fail-Safe Design

Redundant triple valve design

Large Coverage

6" diameter will cover most chest trauma sites

Hydrogel Adhesive

Create a strong seal while allowing repositioning

BCS01

NSN: 6510-01-549-0939

Size: 7.5" x 8" x .25"

Weight: 2.1 oz.

Sold individually, in packs of 10, or cases of 25

Tested to MIL-STD-810G

Shelf Life: 5 years

→ HHMEDCORP.COM

328 McLaws Circle, Williamsburg, VA 23185

(804) 642-3663

hhmedcorp.com



Military Grade. Battle Tested.™

Keeping Someone Alive When Seconds Count.