



H*VENT

Award-Winning Design

- **FAIL-SAFE DESIGN**
6 ports allow for redundancy
- **MULTI-DIRECTIONAL DESIGN**
allow the fluids to “drain with gravity”
- **LARGE PULL TAB**
helps with application and repositioning
- **LARGE VENT OPENINGS**
allow for visual monitoring of the wound

SINGLE PACK (HVS01)
NSN 6510-01-671-7225
Dimensions: 7.75" x 4.5" x 0.25"
Weight: 1.4 oz.

TWIN PACK (HVS02)
Dimensions: 9.25" x 8.5" x 0.25"
Weight: 2.9 oz.



→ HHMEDCORP.COM

328 McLaws Circle, Williamsburg, VA 23185
(804) 642-3663
hhmedcorp.com
2020/07



Military Grade. Battle Tested.™
Keeping Someone Alive When Seconds Count.



Also Available Penetrating Injury Kit

The H*VENT Penetrating Injury Kit combines our newest chest seal with our most popular occlusive chest dressing into a single, easy to open 5 mil pouch for fast access and application.

→ [HHMEDCORP.COM](https://www.hhmedcorp.com)



328 McLaws Circle, Williamsburg, VA 23185
(804) 642-3663
hhmedcorp.com
2020/07



Military Grade. Battle Tested.™
Keeping Someone Alive When Seconds Count.