



TPAK10

TPAK-10

NSN 6515-01-673-1701


Vol. 1000 per case

Total Weight: 15.5kg/34lbs

10 gauge TPAK chest decompression needle is in the recommended 3.25" length for treating a tension pneumothorax. You can't count on a standard IV catheter to needle a chest deeply enough in needle thoracostomy. Our 10 gauge x 3.25" TPAK chest decompression needle is a compact, reliable solution for treating a tension pneumothorax. Shock resistant case ensures the needle/catheter is secure upon deployment, and the capless flash chamber gives instant confirmation of needle placement in pre-hospital needle thoracostomy. TPAK comes in a rugged case, with a textured, easy-open hex cap.

Details:

- 10 gauge x 3.25" chest decompression needle
- Strong, reliable tension pneumothorax needle allows for secure placement of catheter for continuous relief during needle thoracostomy
- Rugged tube with textured, easy-open hex cap
- Capless flash chamber for instant confirmation of needle placement
- Convenient, compact size for maximum portability
- Shock resistant, precision dimensioned case and needle hub ensure that chest decompression needle and catheter secure upon deployment for safety



328 McLaws Circle, Williamsburg, VA 23185

(804) 642-3663

hhmedcorp.com

2020/07



Military Grade. Battle Tested.™

Keeping Someone Alive When Seconds Count.